

Microscope Setup Instructions

1) Unboxing and Checking for Parts

- a) Materials in the manila envelope
 - i) Microscope (with lens cover) and stand
 - ii) MicroUSB adapter
 - iii) Calibration slide

2) Connecting the Camera to Your Computer

- a) *Mac Users*
- b) *Windows Users*
 - i) Open camera app
 - ii) Plug microscope with MicroUSB adapter to MicroUSB port
 - iii) Once connected, the camera should have a flip option, flip the camera
- c) *Other Users*

3) Focusing the Camera

- a) Remove lens cover and angle camera directly upwards
- b) Place blank glass slide on camera and turn knob back and forth until slide is in focus
- c) Once camera is focused replace blank slide with stained slide
- d) Continue to adjust focus knob until cells on slide become visible

4) Powering Down Camera

- a) Turn off camera by unplugging from computer
- b) Be sure to include lens cover when packing up